



## INTRODUCING **Stratus 5**

### The 5 LPM Stationary Oxygen Concentrator from React Health

- Slim line design for a smaller footprint in the home
- 360° swivel wheels to effortlessly roll over all types of flooring, including carpet
- Quiet operation - compressor is specially placed, and padded to reduce noise
- Simple, intuitive operation



REACTHEALTH



**Stratus 5L**  
Designed with the Patient in Mind

# STRATUS 5L

Designed with your business in mind



## Simple to service

Easy access to internal parts

## 3-year standard warranty

### PRODUCT SPECIFICATIONS

|                          |                                |
|--------------------------|--------------------------------|
| Rated power              | 320W                           |
| Voltage                  | 120V $\pm$ 12V, 60HZ $\pm$ 1HZ |
| Flow rate (L/min)        | 0.5 - 5.0                      |
| Oxygen Concentration     | (93% $\pm$ 5%)                 |
| Sound pressure level dBA | ( $\leq$ 45)                   |
| Sound power level dBA    | ( $\leq$ 55)                   |
| Net weight               | 34.1 lbs                       |
| Dimensions               | 15.3 x 9 x 23.6 inches         |



### STRATUS 5 ALARMS

- Abnormal temperature alarm
- Power failure alarm
- Low flow alarm
- Low oxygen purity alarm
- Pressure alarm
- Thermistor fault alarm

### PRODUCT CODES

#### Part Numbers

**STR1005** - Stratus 5 Stationary Oxygen Concentrator

#### Replacement Part Numbers

**STR1020** - Stratus 5 Repl. HEPA Paper Filter- CS of 10

**STR1030** - Stratus 5 Repl. Sponge Filter- CS of 10

#### Oxygen Accessory Part Numbers

**O2S2007** - Standard 7' Oxygen Cannula

**O2U2007** - Ultra Soft 7' Oxygen Cannula

**O2S2012** - Standard 12' Oxygen Cannula

**O2U2012** - Ultra Soft 12' Oxygen Cannula

**O2T2025** - Standard 25' Oxygen Tubing

**O2C2030** - Oxygen Tubing Connector



REACTHEALTH

info@reacthealth.com | www.reacthealth.com