

HomeFill[®]

Oxygen System



The HomeFill Oxygen System by React Health promotes patient ambulation by allowing them to fill their own cylinders from their 5-liter or 10-liter oxygen concentrators.



Provider Benefits

- » Significantly reduces operating costs through fewer cylinder deliveries and service calls and lower storage requirements, all without sacrificing patient care
- » No extra equipment needed because patients can breathe off concentrator while filling cylinders
- » Encourages compliance with oxygen therapy as patients have no worries about running out of oxygen

Features

- » Provides an unlimited, refillable supply of ambulatory oxygen
- » Patient breathes from a continuous flow oxygen concentrator at home
- » Patient protection ensures oxygen purity in the cylinder is always >90%
- » Easy cylinder connections with the OxyLock™ Coupler

Patient Benefits

- » Eliminating oxygen deliveries increases independence, convenience and comfort
- » Small, lightweight cylinders allow for greater mobility than larger, traditional oxygen cylinders
- » Using small, refillable cylinders indoors can help avoid trips and falls associated with using a concentrator with long oxygen tubing

Compatible Concentrators

- » Platinum® 5NXG
- » Platinum 10 Liter
- » Perfecto 2™ Series



Specifications

Dimensions	15" x 20.5" x 16"
Weight	33 lbs.
Electrical Requirements	120VAC +/- 10% @ 60Hz
Limited Warranty	3 Years / 5,000 hours

Part Numbers

HomeFill Compressor	IOH200
Ready Rack for Platinum 10	IOH260
Ready Rack for Perfecto2 Series	IOH270
Ready Rack for Platinum 5NXG	IOH280
M9 Cylinder w/ Integrated Conserver	HF2PC9
M9 Cylinder w/ Post Valve	HF2POST9
D Cylinder w/ Post Valve	HF2POSTD